

Disinfection Cabinet

DC-430



Easy Banker International Co., Ltd
www.easy-banker.com

Product use and characteristics

The banknote disinfection cabinet adopts a combination of ultraviolet and ozone disinfection factor disinfection methods. Through the ultraviolet radiation and the strong oxidative decomposition ability of ozone, it destroys the internal structure of bacteria and virus DNA, making it immediately die and the market reproduction ability. For cash collected at bank outlets and cash deposited in financial self-service equipment, and cash collected by financial institutions from hospitals, farmers' markets, buses, enterprises, etc., use the cash disinfection cabinet to complete the sterilization work,

Disinfection method

ozone + ultraviolet

Specifications

Model:	DC-430
Cabinet volume:	430L
Dimensions:	56*47*170cm
Packaging parameters:	70*62*190cm
Power supply:	220V-50HZ
Rated input power:	150W-1800W
Net weight:	35kg
Gross weight:	37kg
Cash capacity(sheets):	50000
Screen load:	20KG
Service life≥1000h:	UV generation life≥1000h
Material:	Brushed stainless steel
DC-50Material:	Paint+Tempered glass
Full load index parameters:	Ozone concentration Disinfection time 40mg/m ³ > 0m ³ ≥20mg/m ³ 0-60min



Easy Banker International Co., Ltd
www.easy-banker.com