

Disinfection Cabinet

DC-300



Easy Banker International Co., Ltd
www.easy-banker.com

Product use and characteristics

The banknote disinfection cabinet adopts a combination of ultraviolet and ozone disinfection factor disinfection methods. Through the ultraviolet radiation and the strong oxidative decomposition ability of ozone, it destroys the internal structure of bacteria and virus DNA, making it immediately die and the market reproduction ability. For cash collected at bank outlets and cash deposited in financial self-service equipment, and cash collected by financial institutions from hospitals, farmers' markets, buses, enterprises, etc., use the cash disinfection cabinet to complete the sterilization work,

Disinfection method

ozone + ultraviolet

Specifications

Model:	DC-300
Cabinet volume:	300L
Dimensions:	56*47*125cm
Packaging parameters:	70*62*144cm
Power supply:	220V-50HZ
Ratedinputpower:	150W-1800W
Net weight:	30kg
Gross weight:	32kg
Cash capacity(sheets):	40000
Screen load:	20KG
Service life \geq 1000h:	UV generation life \geq 1000h
Material:	Brushed stainless steel
DC-50Material:	Paint+Tempered glass
Full load index parameters:	Ozone concentration Disinfection time 40mg/m ³ > 0m ³ \geq 20mg/m ³ 0-60min



Easy Banker International Co., Ltd
www.easy-banker.com